

IN APPLICATION

Progress in Volumetric Flow Measurement

LaVision's FlowMaster Tomo PIV applied to a turbulent jet in water at the TU Delft, Aerospace Engineering, Aerodynamic Section



Introduction

Conventional stereo PIV only allows the measurement of cross sections in three dimensional flow fields. Volumetric flow measurement enables time resolved 3-dimensional, 3-component velocity measurement of complete flow structures.

Recent progress of the volumetric flow measurement system FlowMaster Tomo PIV from LaVision is demonstrated using experimental data recorded at the TU Delft.

The progress includes three main features that are exemplified in this application note:

1. **Fast preview mode**
2. **3 camera Tomo PIV**
3. **Motion Tracking Enhanced reconstruction**


Experimental Setup

Time resolved jet flow in water at $Re = 5000$ is studied. The jet diameter is 10 mm. A high speed laser (2×25 mJ/pulse) illuminates a cylinder of 50 mm diameter. Four HighSpeedStar 6 cameras ($1k \times 1k$) record instantaneous volumes of $50 \times 50 \times 60$ mm³ at 1.3 kHz recording rate.



picture from jet setup with LaVision's new Imager pro HS 4M cameras

Recordings at different seeding particle densities are used in the following to demonstrate the fast preview mode (low density: 0.031 particles per pixel), 3 camera Tomo PIV (moderate density: 0.055 particles per pixel) and motion tracking enhanced reconstruction (highest density: 0.17 particles per pixel).



LaVisionUK Ltd

2 Minton Place / Victoria Road
Bicester, Oxon / OX26 6QB / United Kingdom
E-Mail: sales@lavision.com / www.lavisionuk.com
Phone: +44-(0)-870-997-6532 / Fax: +44-(0)-870-762-6252

LaVision GmbH

Anna-Vandenhoeck-Ring 19
D-37081 Göttingen / Germany
E-Mail: info@lavision.com / www.lavision.com
Tel. +49-(0)551-9004-0 / Fax +49-(0)551-9004-100

LaVision Inc.

211 W. Michigan Ave. / Suite 100
Ypsilanti, MI 48197 / USA
E-mail: sales@lavisioninc.com / www.lavisioninc.com
Phone: (734) 485 - 0913 / Fax: (240) 465 - 4306



Progress in Volumetric Flow Measurement

LaVision's FlowMaster Tomo PIV applied to a turbulent jet in water at the TU Delft, Aerospace Engineering, Aerodynamic Section

Fast preview mode

Low particle density: 0.031 particles per pixel, 4 cameras.

The fast **Multiplicative Line of Sight reconstruction (MLOS)**, together with only 50% overlap for correlation (interrogation volume size 48^3 volume elements = voxels) allows a rapid preview of the flow structure. $25 \times 40 \times 25$ vectors already show the principal flow features after a computation time of only 85 seconds. [MLOS reconstruction ($891 \times 965 \times 595$ voxels) + correlation (2 passes)].



3 camera Tomo PIV

Moderate particle density: 0.055 particles per pixel

Volume Self Calibration¹⁾ reduces the ghost intensity at a given particle density. This allows the usage of only three cameras at moderate seeding densities.


Standard **Multiplicative Algebraic Reconstruction**

Technique (MART)²⁾ is applied for the computation of well resolved flow fields at moderate computation times. The interrogation volume size is again 48^3 voxels. 75% overlap results in $51 \times 81 \times 51$ vectors.

Motion Tracking Enhanced reconstruction

High particle density: 0.17 particles per pixel, 4 cameras

Motion Tracking Enhanced MART (MTE)³⁾ is the reconstruction technique that allows the highest seeding particle density. Due to the high density of 0.17 particles per pixel, the interrogation volume size can be reduced to 32^3 voxels. This reduction of the interrogation volume size from 48^3 to 32^3 results in $1.5^3 = 3.4$ times more vectors: $76 \times 121 \times 76$, increasing the spatial resolution so that even small flow structures can be examined in more detail. **MTE** also leads to a significant increase of the velocity field accuracy⁴⁾.



¹⁾Wieneke B. (2008) Volume Self-Calibration for 3D Particle Image Velocimetry. Exp in Fluids 45: 549-556

²⁾Elsinga G. E., Scarano F., Wieneke B., van Oudheusden B.W. (2006) Tomographic particle image Velocimetry. Exp in Fluids 41: 933-947

³⁾Novara M., Batenburg K.J., Scarano F. (2010) Motion tracking-enhanced MART for tomographic PIV. Meas. Sci. Technol. 21 035401

⁴⁾Michaelis D., Novara M., Scarano F., Wieneke B. (2010) Comparison of volume reconstruction techniques at different particle densities, 15th Int Symp on Applications of Laser Techniques to Fluid Mechanics, Lisbon, Portugal