

Helium-filled Soap Bubble Generator

for large volume PIV & Particle Tracking Helium-filled Soap Bubbles (HFSB) were successfully used in the past for flow visualization and large scale PIV & PTV applications in air [1]. The mm-sized bubbles are neutrally buoyant and extremely bright. The main limitation was the low production rate, e.g. <1000 per sec., not sufficient for many applications.



LaVision now developed in cooperation with DLR Göttingen and TU Delft a new seeding generator and nozzle design for up to 200 nozzles in parallel. The bubbles have a highly constant size of down to 300 µm and are 10,000 times brighter than standard oil aerosol seeding particles. With a

Bubbles leaving nozzle at a rate of 40,000/sec.

relaxation time of 11 µs they follow the flow perfectly, even at high velocities and flow gradients [2]. Large measurement volumes can be filled with a high seeding concentration, recorded and processed with LaVision's time-resolved Tomographic PIV or Shake-the-Box particle tracking systems.

Micro-processor controlled operation

Multi-nozzle arrays

The micro-processor inside the controller device automatically manages the air, helium and soap flow rates with start, stop and pause modes. All maintenance procedures such as cleaning and filling are carried out at a push button. A large soap reservoir enables continuous operation of up to 20 h. Remote control is possible with LAN connection to any remote device with standard internet browser.

The nozzles are separated from the controller at a distance of up to 20 m. An aerodynamically shaped Linear Nozzle Array (LNA) with up to 20 nozzles has been designed for e.g. wind tunnel operation, typically placed inside the settling chamber. Several LNAs can be connected to a single controller for flexible configurations and seeding large measurement volumes of more than 1 m³.



LaVisionUK Ltd

2 Minton Place / Victoria Road Bicester, Oxon / OX26 6QB / United Kingdom E-Mail: sales@lavision.com / www.lavisionuk.com Phone: +44-(0)-870-997-6532 / Fax: +44-(0)-870-762-6252

LaVision GmbH


LaVision Inc.

Anna-Vandenhoeck-Ring 19 D-37081 Göttingen / Germany E-Mail: info@lavision.com / www.lavision.com Tel. +49-(0)551-9004-0 / Fax +49-(0)551-9004-100

211 W. Michigan Ave. / Suite 100 Ypsilanti, MI 48197 / USA E-mail: sales@lavisioninc.com / www.lavisioninc.com Phone: (734) 485 - 0913 / Fax: (240) 465 - 4306




Example



Volume: 800 x 800 x 1200 mm³, illuminated with LED spots

Thermal air convection in 1 cubic meter



Particle tracks and vorticity contours calculated with Shake-the-Box (STB), courtesy D. Schanz, A. Schröder, F. Huhn DLR Göttingen

Specifications

bubble size: down to 300 μm

- 10,000 brighter compared to usual oil aerosol seeding
- LED illumination possible in many applications
- neutrally buoyant, relaxation time $\tau = 11 \mu s$
- bubble lifetime of several minutes
- > aerodynamically shaped Linear Nozzle Array (LNA) with up to 20 nozzles
- large measurement volume > 1 m³
- high seeding concentration: 1300 bubbles per cm³ at nozzle exit
- > remote control via LAN: start, stop, pause, etc.

References:

[1] Bosbach J, Kühn M, Wagner C (2009) Large scale particle image velocimetry with helium-filled soap bubbles. Exp Fluids 46:539–547

[2] Scarano F, Sina Ghaemi S, Carlo Alp Caridi G, Bosbach J, Dierksheide U, Sciacchitano A (2015), On the use of helium-filled soap bubbles for large-scale tomographic PIV in wind tunnel experiments, Exp Fluids 56:42

Data provided by LaVision are believed to be true. However, no responsibility is assumed for possible inaccuracies or omissions. All data are subject to change without notice.

Dec-18

LaVisionUK Ltd

2 Minton Place / Victoria Road Bicester, Oxon / OX26 608 / United Kingdom E-Mail: sales@lavision.com / www.lavisionuk.com Phone: +44-(0)-870-997-6532 / Fax: +44-(0)-870-762-6252

LaVision GmbH

LaVision Inc.

Anna-Vandenhoeck-Ring 19 D-37081 Göttingen / Germany E-Mail: info@lavision.com / www.lavision.com Tel. +49-(0)551-9004-0 / Fax +49-(0)551-9004-100 211 W. Michigan Ave. / Suite 100 Ypsilanti, MI 48197 / USA E-mail: sales@lavisioninc.com / www.lavisioninc.com Phone: (734) 485 - 0913 / Fax: (240) 465 - 4306